

Bst DNA Polymerase (Large Fragment)

v. 230101

Catalog number C15019-1600U / C15019-8000U

Package & Component	Cat.	Name	Amount	
	C15019-1600U		Bst DNA Polymerase (Large Fragment) (8 U/ μ L)	1,600 U
			10X Bst DNA Polymerase Reaction Buffer	1 mL
			100 mM MgSO ₄	0.4 mL
	C15019-8000U		Bst DNA Polymerase (Large Fragment) (8 U/ μ L)	8,000 U
			10X Bst DNA Polymerase Reaction Buffer	1 mL
		100 mM MgSO ₄	0.4 mL	

Description Bst DNA Polymerase (Large fragment) is an enzyme of *Bacillus stearothermophilus* DNA polymerase which can catalyze 5' \rightarrow 3' polymerase activity but lacks 5' \rightarrow 3' exonuclease activity. Bst DNA Polymerase offers strand displacement capabilities, making it ideal for isothermal amplification.

Source *Escherichia coli*

Purity >98% as determined by SDS-PAGE (purified by Ni-NTA chromatography).

Unit Definition One unit is defined as the amount of the enzyme incorporates 10 nmol of dNTP into acid-insoluble product in 30 minutes at 65°C.

Reaction Condition 1X Bst DNA Polymerase reaction buffer, supplemented with dNTP mix, primer and DNA template. Incubate at 65°C.
10X Bst DNA Polymerase Reaction Buffer: 200 mM Tris-HCl (pH 8.8), 100 mM (NH₄)₂SO₄, 100 mM KCl, 20 mM MgSO₄, and 1% Triton® X-100.

Storage Buffer Bst DNA Polymerase (Large fragment) is supplied in 10 mM Tris-HCl (pH7.5), 50 mM KCl, 1 mM DTT, 0.1 mM EDTA, 0.1% Triton® X-100 and 50% Glycerol.

Storage Stored at -20°C. Avoid repeated freeze/thaw cycles.

LAMP reaction recipe:

1. Place all required reagents **on ice** and add each of them following the order suggested below.

Component	Amount	Final concentration
10X Bst DNA Polymerase Reaction Buffer	2.5 μ L	1X
100 mM MgSO₄	1.5 μ L	6 mM final concentration, total 8 mM
10 mM dNTP mix	3.5 μ L	1.4 mM each
10X FIP/BIP primers	1 μ L	1.6 μ M
10X F3/B3 primers	1 μ L	0.2 μ M
10X LoopF/B primers	1 μ L	0.8 μ M
DNA template	X μ L	10 copies or more
Nuclease-Free H₂O	Y μ L	-
Bst DNA Polymerase (Large Fraction) (8 U/μL)	1 μ L	8 U/rxn
Total reaction volume	25 μ L	-

Manuel

2. Gently mix the reaction thoroughly to achieve uniform distribution.
3. Incubate at 65°C for 30-60 minutes.
4. MgSO₄ (2-10 mM), Bst DNA Polymerase (40-320 U/mL) and temperature (50-65 °C) can be adjusted for optimal results.
5. Reaction preparations may be scaled up or down proportionately.

Notes

It is not recommended to perform reaction above 70 °C. Bst DNA Polymerase cannot be used for thermal cycle sequencing.

For Research Use Only.