

**HIV-1 RTase (HIV-1 Reverse Transcriptase)**

v. 230101

<b>Catalog number</b>	C15020-200U		
<b>Package &amp; Component</b>	<b>Cat.</b>	<b>Name</b>	<b>Amount</b>
	C15020-200U	HIV-1 Reverse Transcriptase (5 U/μL)	200 U
		10X Isothermal Amplification Buffer	1 mL
		100 mM MgSO <sub>4</sub>	0.4 mL
<b>Description</b>	<p>The human immunodeficiency virus type 1 reverse transcriptase (HIV-1 RTase is an enzyme that can catalyze complementary DNA (cDNA) synthesis from an RNA template. Due to its greater thermostability than comparatives of AMV and MMLV, HIV-1 RTase is currently used for RT-LAMP reactions, in combination with Bst DNA polymerase LF.</p>		
<b>Source</b>	Escherichia coli		
<b>Purity</b>	>98% as determined by SDS-PAGE (purified by Ni-NTA chromatography).		
<b>Unit Definition</b>	One unit is defined as the amount of the enzyme incorporates 1 nmol of dTTP into acid-insoluble product in 10 minutes at 50°C.		
<b>Reaction Condition</b>	<p>1X Isothermal Amplification Buffer, supplemented with dNTP mix, template and primer, Incubate at 65°C for one-step RT-LAMP.          10X Isothermal Amplification Buffer: 200 mM Tris-HCl (pH 8.8), 100 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>, 500 mM KCl, 20 mM MgSO<sub>4</sub>, and 1% Tween 20.</p>		
<b>Storage Buffer</b>	HIV-1 Reverse Transcriptase is supplied in 10 mM Tris-HCl (pH 7.4), 100 mM KCl, 1 mM DTT, 0.1 mM EDTA and 50% glycerol (v/v).		
<b>Storage</b>	Stored at -20°C. Avoid repeated freeze/thaw cycles.		

**First strand cDNA synthesis:**

- Place all required reagents **on ice** and add each of them following the order suggested below.

	<b>Component</b>	<b>Amount</b>	<b>Final concentration</b>
<b>Manuel</b>	<b>Oligo (dT)12-18 (50 μM) or random primer mix (60 μM)</b>	1 μL	-
	<b>Total RNA template</b>	X μL	1 μg
	<b>Nuclease-Free H<sub>2</sub>O</b>	Y μL	-
	<b>Total reaction volume</b>	13.5 μL	-

- Heat the tube to 65°C for 10 minutes to denature the secondary structure within RNA template. Immediately cool the tube on ice for 1 minute and centrifuge briefly in microcentrifuge. Add the following components to the annealed primer/RNA template, prepare on ice.

Component	Amount	Final concentration
<b>5X Reverse Transcriptase Reaction Buffer</b>	4 µL	1X
<b>10 mM dNTPs (each)</b>	1 µL	0.5 mM each
<b>RNase Inhibitor</b>	0.5 µL	20 U/rxn
<b>HIV RTase</b>	1 µL	5 U/rxn
<b>Total reaction volume</b>	µL	-

- Incubate at 42-50°C for 1 hour.
- Inactivate the reaction at 80°C for 10 minutes. The cDNA products should be store at -20°C.
- Reaction preparations may be scaled up or down proportionately.

#### **RT-LAMP reaction.**

- Place all required reagents **on ice** and add each of them following the order suggested below.

Component	Amount	Final concentration
<b>10X Isothermal Amplification Buffer</b>	2.5 µL	1X
<b>100 mM MgSO<sub>4</sub></b>	1.5 µL	6 mM final concentration, total 8 mM
<b>10 mM dNTPs mix</b>	3.5 µL	1.4 mM each
<b>10X FIP/BIP primers</b>	1 µL	1.6 µM
<b>10X F3/B3 primers</b>	1 µL	0.2 µM
<b>10X LoopF/B primers</b>	1 µL	0.8 µM
<b>RNA template</b>	X µL	-
<b>Nuclease-Free H<sub>2</sub>O</b>	Y µL	-
<b>Bst DNA Polymerase (Large Fraction)</b>	8 U	8 U/rxn
<b>HIV-1 Reverse Transcriptase</b>	2 µL	10 U/rxn
<b>Total reaction volume</b>	25 µL	-

- Gently mix the reaction thoroughly to achieve uniform distribution.
- Incubate at 65°C for 30-60 minutes.
- MgSO<sub>4</sub> (2-10 mM), Bst DNA Polymerase (40-320 U/mL), HIV-1 RT and temperature (50-65 °C) can be adjusted for optimal results.

5. Reaction preparations may be scaled up or down proportionately.

---

**Notes**

After the reaction is complete, Bst DNA Polymerase and HIV-1 RT can be inactivated by incubation at 80°C for 10 minutes.

---

*For Research Use Only.*