

M-MLV RTase (M-MLV Reverse Transcriptase)

v. 230101

Catalog number	C15021-20000U / C15021-50000U		
Package & Component	Cat.	Name	Amount
	C15021-20000U	M-MLV Reverse Transcriptase (200 U/ μ L)	20,000 U
		5X M-MLV Reverse Transcriptase Reaction Buffer	1 mL
	C15021-50000U	M-MLV Reverse Transcriptase (200 U/ μ L)	50,000 U
5X M-MLV Reverse Transcriptase Reaction Buffer		1 mL	
Description	Moloney Murine Leukemia Virus (M-MLV) Reverse Transcriptase is an RNA-dependent DNA polymerase that synthesizes the first strand of complementary cDNA from a single-stranded RNA template with hybridized primer. This kit features high activity formulation of M-MLV RT and 5X reverse transcription buffer which are capable of full-length cDNA synthesis and high cDNA yields.		
Source	Escherichia coli		
Purity	>98% as determined by SDS-PAGE (purified by Ni-NTA chromatography).		
Unit Definition	One unit is defined as the amount of the enzyme incorporates 1 nmol of dTTP into acid-insoluble product in 10 minutes at 37°C.		
Reaction Condition	1X M-MLV Reverse Transcriptase Buffer, supplemented with dNTP mix, template and primer, Incubate at 37°C for synthesis of first strand cDNA. 5X M-MLV Reverse Transcriptase Buffer: 250 mM Tris-HCl (pH 8.3), 15 mM MgCl ₂ , 375 mM KCl, and 50 mM DTT.		
Storage Buffer	M-MLV Reverse Transcriptase is supplied in 20 mM Tris-HCl (pH 7.5), 200 mM NaCl, 0.25 mM EDTA, 0.01% NP-40 (v/v), 2.5 mM DTT, 50% glycerol (v/v).		
Storage	Stored at -20°C. Avoid repeated freeze/thaw cycles.		

First strand cDNA synthesis:

- Place all required reagents to a nuclease-free microcentrifuge tube and following the order suggested below.

Component	Amount	Final concentration
Oligo (dT)12-18 (50 μM) or random primer mix (60 μM)	1 μ L	-
Total RNA template	X μ L	1 μ g
Nuclease-Free H₂O	Y μ L	-
Total reaction volume	10 μ L	-

- Heat the tube to 65°C for 10 minutes to denature the secondary structure within RNA template. Immediately cool the tube on ice for 1 minute and centrifuge briefly in microcentrifuge. Add the following components to the annealed primer/RNA template, prepare on ice.

Manuel

Component	Amount	Final concentration
5X Reverse Transcriptase Reaction Buffer	4 μ L	1X
10 mM dNTPs mix	1 μ L	0.5 mM each
RNase Inhibitor	X μ L	20 U/rxn
M-MLV RTase	1 μ L	200 U/rxn
Nuclease-Free H₂O	Y μ L	-
Total reaction volume	20 μ L	-

- Incubate at 37°C for 1 hour. The extension temperature may be adjusted from 37°C to 42°C.
- Inactivate the reaction at 65°C for 20 minutes. The cDNA products should be store at -20°C.
- Reaction preparations may be scaled up or down proportionately.

Notes

After the reaction is complete, M-MLV RTase can be inactivated by incubation at 65°C for 20 minutes.

For Research Use Only.