

Catalog number	C15028-1ML		
Set package	Cat.	Name	Amount
	C15028-1ML	StablePlus™ 2X RT-LAMP Master Mix	1 mL
Description	<p>Croyez StablePlus™ 2X RT-LAMP Master Mix is an optimized master mix for reverse-transcription loop-mediated isothermal amplification (RT-LAMP) reactions. This product is a dual enzyme system, providing a rapid and sensitive detection in one pot. The amplified products can be detected by agarose gel electrophoresis.</p> <p>The StablePlus™ version contains nucleic acid stabilizing agent to protect the amplified products.</p>		
Storage	Stored at -20°C. Avoid repeated freeze/thaw cycles.		

The following procedure is a general guideline for RT-LAMP reaction. To maintain an RNase-free environment, always wear disposable gloves, and use laboratory consumables and water of nuclease-free grade during the whole experiment course. **RT-LAMP reaction set-up:**

- 10X LAMP primer mix

Component	10X concentration	Final concentration
FIP	16 µM	1.6 µM
BIP	16 µM	1.6 µM
F3	2 µM	0.2 µM
B3	2 µM	0.2 µM
LOOP F	8 µM	0.8 µM
LOOP B	8 µM	0.8µM

Manuel

- An overview of the reaction set-up is listed in the table below. Place all required reagents **on ice**. Distribute appropriate volumes into each tube before adding template.

Component	Amount	Final concentration
StablePlus™ 2X RT-LAMP Master Mix	12.5 µL	1X
10X LAMP primer mix	2.5 µL	1X
RNA template	1-2 µL	variable
Nuclease-Free H₂O	X µL	-
Total reaction volume	µL	-

* See Usage Notes for additional guidelines on primer/template preparation.

3. Add target RNA template to the detection tube. Gently mix the reaction thoroughly to achieve uniform distribution and avoid making bubbles.
4. Incubate at 65°C for 30-60 min.
5. After LAMP reaction complete, the enzyme can be inactivated by heating at 80°C for 10 min.

Usage Notes***Primer concentration***

Primer concentration can be titrated between 0.25X – 1X if undesired background signal appeared.

Detection method

Both detecting the amplified products by agarose gel electrophoresis and turbidity changes due to magnesium pyrophosphate precipitation can be employed to analyze test results, but the latter is somehow less sensitive.

For Research Use Only.