


PRODUCT INFORMATION

**UniHRV 3C Protease**

v. 230201

<b>Catalog number</b>	C09003-bulk / C09003-500U UniHRV 3C Protease C09002-B-bulk / C09002-B-1 10X HRV 3C Cleavage Buffer		
<b>Package</b>	Cat.	Name	Amount
	C09003-500U	UniHRV 3C Protease	500U
		Cleavage Validated Protein	20 µg
		10X HRV 3C Cleavage Buffer	10 mL
<b>Description</b>	UniHRV 3C Protease is a recombinant form of the 3C protease derived from human rhinovirus 14 expressed in <i>E. coli</i> (specific activity 1800-2000 U/mg). This product is a highly purified recombinant 6XHis-fusion protein and requires neither metal nor cofactors for activity. UniHRV 3C Protease recognizes the cleavage site: Leu-Glu-Val-Leu-Phe-Gln↓Gly-Pro (LEVLFG↓GP). UniHRV 3C Protease demonstrate excellent cleavage efficiency in a variety of fusion proteins.		
<b>Source</b>	<i>Escherichia coli</i>		
<b>Endotoxin level</b>	<1 EU per 1 µg of the protein by the LAL method.		
<b>Activity</b>	One unit of HRV 3C Protease is defined as the amount of enzyme needed to digest 1 nM of the pNA-peptide substrate per 10 min at RT or per hour at 0 °C in the reaction buffer.		
<b>Purity</b>	>98% as determined by SDS-PAGE analysis.		
<b>Formulation</b>	UniHRV 3C Protease was lyophilized from a solution containing 50 mM Tris, 150 mM NaCl, 1 mM EDTA, 1 mM DTT, 0.04% Tween20, 8% trehalose, 8% mannitol. 10X HRV 3C Cleavage Buffer: 1.5 M NaCl, 0.5 M Tris-HCl, pH 7.5		
<b>Reconstitution</b>	It is recommended to reconstitute the lyophilized protein in sterile H <sub>2</sub> O and incubate the stock solution for at least 20 min to ensure sufficient re-dissolved.		
<b>Storage</b>	Lyophilized protein should be stored at -20°C. Upon reconstitution, protein aliquots should be stored at -20°C. HRV 3C Protease Cleavage Buffer should be stored at -20°C or 4°C.		



### SDS-PAGE analysis of UniHRV 3C Protease

The following protocol is an example of a simple optimization experiment designed to estimate the appropriate ratio of enzyme:target protein. This example represents UniHRV 3C Protease:target protein ratios (U/ $\mu$ g) of 1:10, 1:25 and 1:50.

Component	Volume ( $\mu$ L)
UniHRV 3C Protease (1 U/ $\mu$ L)	1
Validated Protein (10, 25, and 50 $\mu$ g each)	X
10X HRV 3C Cleavage Buffer	5
H <sub>2</sub> O	44-X
Total volume	50

- Incubate the reaction mixture at 4°C for 16 hours or overnight.
- Determine cleavage level of the samples by SDS-PAGE analysis.
- UniHRV 3C protease: target protein ratio of 1:5~1: 200 (U/ $\mu$ g) is used for most fusion protein cleavage.
- If shorter incubation time is required, more amount of UniHRV 3C protease or higher temperature (RT) can be implemented. Reaction can be performed at 4°C-37°C. 4°C is recommended as the starting standard.

### Manuel




Fig. The validated protein was cleaved by UniHRV 3C protease and Competitor's HRV 3C Protase at 4°C for 16 h.

---

**Notes (Important)**

- Cleavage efficiency may differ based on structure and properties of each target protein, we recommend testing several enzyme-to-substrate ratios, temperatures, and incubation times.
- HRV 3C Protease reactions can be performed in a buffer which is optimal for the target protein. Reducing reagents (e.g., DTT) or salts (e.g., NaCl) can be added for cleavage efficiency evaluation.
- Usually, 1 U UniHRV 3C Protease can cleavage >95% of target protein (from 50 to 400 µg) at 4°C for 16h. Even different lot of target proteins might result in different cleavage efficiency. The validated protein (50 µg) provide in this product can be stably cleavage >50% when using 1 U of UniHRV 3C Protease.

---

*For Research Use Only. Not for use in diagnostic and therapeutic procedures.*