

PRODUCT INFORMATION

T7 RNA Polymerase Transcription Buffer Set

v. 230501

Catalog number C15027-K01 / C15027-K02

Package & Component		C15027-K01	C15027-K02
	T7 RNA Polymerase (200 U/μL)	10,000 U	25,000 U
	10X RNA Polymerase Reaction Buffer A	0.5 mL	1 mL
	10X RNA Polymerase Reaction Buffer B	0.5 mL	1 mL
	10X RNA Polymerase Reaction Buffer C	0.5 mL	1 mL
	10X RNA Polymerase Reaction Buffer D	0.5 mL	1 mL
	10X RNA Polymerase Reaction Buffer E	0.5 mL	1 mL
	10X RNA Polymerase Reaction Buffer F	0.5 mL	1 mL
	10X RNA Polymerase Reaction Buffer G	0.5 mL	1 mL
	100 mM DTT	0.5 mL	1 mL

Description

The T7 RNA Polymerase Transcription Buffer Set contains a T7 RNA Polymerase and a set of seven 10X RNAP reaction buffers. This set of reagents was specifically designed for the selection of efficient transcription reactions. The user simply combines DNA template, NTPs, T7 RNA polymerase, and the 10X polymerase reaction buffer to proceed following reactions.

Standard RNA synthesis procedures:

1. Below reaction mixture should be prepared under room temperature and combined in the following order:

Manuel

Component	Amount	Final concentration
Nuclease-Free H₂O	X μL	-
Template DNA	0.5-1 μg	-
10X RNA Polymerase Reaction Buffer	2 μL	1 X
ATP (100 mM)	1 μL	5 mM
UTP (100 mM)	1 μL	5 mM
CTP (100 mM)	1 μL	5 mM
GTP (100 mM)	1 μL	5 mM
100 mM DTT	2 μL	10 mM
T7 RNA Polymerase (200 U/μL)	1 μL	-
RNase inhibitor (optional)	0.5 μL	1 U/μL
Total reaction volume	20 μL	-

-
2. Incubate at 37°C for 30 minutes to 2 hours.
 3. Above reaction mixture may be scaled up or down proportionately.
-

Storage	Stored at -20°C. <i>For optimal storage, aliquot the enzyme, reaction buffer and DTT reagent into smaller quantities and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.</i>
----------------	--

- | | |
|--------------|--|
| Notes | <ol style="list-style-type: none">1. Transcription reaction should be performed under RNase free condition. Use nuclease-free tubes, reagents, and water to avoid RNase contamination. Also, wear gloves when working with RNA.2. To obtain optimal condition, NTP concentration can be titrated between 5 – 10 mM.3. The volume of T7 RNA Polymerase can be titrated between 1-2 µL in the IVT reaction to optimize your assay. |
|--------------|--|
-

For Research Use Only.