

PRODUCT INFORMATION

GDNF (Glial-derived neurotrophic factor), Mouse

Catalog number	C02089-5UG / C02089-20UG / C02089-100UG
Package	5 µg / 20 µg / 100 µg
Description	Glial cell-derived neurotrophic factor (GDNF) is a protein that, in humans, is encoded by the GDNF gene. GDNF is a small protein that potently promotes the survival of many types of neurons. GDNF, that acts via classical neurotrophic mechanism, has been effective in several pre-clinical models of PD and had some efficacy in parkinsonian patients.
Source	<i>Escherichia coli</i>
Sequence	MSPDKQAAALPRRERNRQAAAASPENSRGKGRRGQRGKNRGCVLTAIHLNVT DLGLGYETKEELIFRYCSGSCESAETMYDKILKNLSRSRRLTSDKVGQACCRPV AFDDDL SFLDDNLVYHILRKHS AKRCGCI with polyhistidine tag at the C-terminus
Endotoxin level	<0.1 EU per 1 µg of the protein by the LAL method.
Purity	>98% as determined by SDS-PAGE. Ni-NTA chromatography
Formulation	The protein was lyophilized from a solution containing 1X PBS, pH 7.4.
Reconstitution	It is recommended to reconstitute the lyophilized protein in sterile H ₂ O to a concentration not less than 100 µg/mL and incubate the stock solution for at least 20 min to ensure sufficient re-dissolved.
Storage	Lyophilized protein should be stored at -20°C. Upon reconstitution, protein aliquots should be stored at -20°C or -80°C.
Note	Please use within one month after protein reconstitution.



SDS-PAGE analysis of recombinant mouse GDNF

For Research Use Only.