

PRODUCT INFORMATION

**SARS-CoV-2 IgG Neutralizing Antibody Rapid Test (Version II) /
UNCUT SHEET**

Catalog number	C00003
Package	Sheet
Description	SARS-CoV-2 IgG Neutralizing Antibody Rapid Test provides a simple, convenient testing for SARS-CoV-2 IgG neutralizing antibody of serum. The neutralizing antibody blocks the interaction between RBD and ACE2. The uncut sheet must be further cut into individual strips. (Test line: anti-human IgG)
Product type	Uncut sheet
Reactivity	SARS-CoV-2
Storage	Store at 15-30°C and keep dry for one year Avoid direct exposure to sunlight Do not freeze
Application	LFIA

For Research Use Only.

