

PRODUCT INFORMATION
GTP Solution (100 mM)

Catalog number	C15025-1ML		
Set package	Cat.	Name	Amount
	C15025-1ML	GTP (100 mM)	1 mL
Description	Croyez Guanosine-5'-triphosphate (GTP) is supplied as a 100 mM aqueous solution at pH 7.5.		
General characteristics	GTP (Guanosine-5'-triphosphate, Sodium salt)		
Purity	≥ 98% purity of each component confirmed by HPLC.		
Form	Clear aqueous solutions (pH 7.4).		
Storage	Stored at -20°C. Avoid repeated freeze/thaw cycles.		
Application	In vitro transcription		
	RNA amplification		
	miRNA and siRNA synthesis		
Notes	Mix thoroughly before use and dilute with nuclease-free water if necessary.		

For Research Use Only.

