

## PRODUCT INFORMATION

**Goat anti-Mouse IgG (H+L)-HRP**

v. 230201

<b>Catalog number</b>	C04001-2ML / C04001-100ML
<b>Package</b>	2x1 mL / 100 mL
<b>Description</b>	Croyez Goat anti-Mouse IgG (H+L)-HRP is conjugated with high-activity HRP. In direct ELISA, less than 1% cross-reactivity with human, goat, and rabbit IgG is observed.
<b>Product type</b>	Polyclonal
<b>Concentration</b>	Lot-dependent, please refer to product label.
<b>Immunogen</b>	Purified mouse IgG
<b>Host</b>	Goat
<b>Reactivity</b>	Mouse (Minor cross-reactivity to rat)
<b>Conjugation</b>	HRP (RZ>3)
<b>Isotype</b>	N/A
<b>Purity</b>	Immunogen affinity chromatography.
<b>Form</b>	Liquid
<b>Storage buffer</b>	Phosphate Buffered Saline containing 1% BSA, 0.01% thimerosal, and 50% glycerol, pH 7.4.
<b>Storage</b>	Store at 4°C for at least 6 months. Store at -20°C for long-term preservation. Suggestion: Divide antibody into several vials. Keep only vials for usage at 4°C.
<b>Application</b>	ELISA, IHC, WB
<b>Application Note</b>	N/A
<b>Manual</b>	1:500-1:2500 for IHC 1:5000-1:25000 for ELISA and WB Working dilution for specific application should be determined by the investigator to obtain the best conditions.

*For Research Use Only.*