


PRODUCT INFORMATION

SARS-CoV-2 Spike Protein RBD, His-SUMO Tag, HEK293

v. 230201

Catalog number	C11020-100UG
Package	100 µg
Source	HEK293
Sequence	RVQPTESIVRFPNITNLCPFGEVFNATRFASVYAWNRRKRISNCVADYSVLYNSA SFSTFKCYGVSPTKLNLDLCTNVYADSFVIRGDEVQRQIAPGQTGKIADYNYKLP DDFTGCVIAWNSNNLDSKVGGNYNLYRLFRKSNLKPFFERDISTEIQAGSTPC NGVEGFNCYFPLQSYGFQPTNGVGYPYRVVVLSFELLHAPATVCGPKKSTNL VKNKCVNF with polyhistidine-SUMO tag at the N-terminus
Endotoxin level	<0.1 EU per 1 µg of the protein by the LAL method.
Purity	>98% as determined by SDS-PAGE analysis. Purified by Ni-NTA chromatography.
Form	Liquid
Storage buffer	Phosphate buffered saline, pH 7.4.
Storage	Stored at -80°C. This product should be used and stored under sterile condition. Avoid repeated free-thaw cycles.



SDS-PAGE analysis of SARS-CoV-2 Spike Protein RBD, His-SUMO Tag, HEK293

For Research Use Only.