

PRODUCT INFORMATION

ATP Solution (100 mM)

v. 230101

Catalog number	C15022-1ML		
Set package	Cat.	Name	Amount
	C15022-1ML	ATP (100 mM)	1 mL
Description	Croyez Adenosine 5'-triphosphate (ATP) is supplied as a 100 mM aqueous solution at pH 7.5.		
General characteristics	ATP (Adenosine-5'-triphosphate, Sodium salt)		
Purity	≥ 98% purity of each component confirmed by HPLC.		
Form	Clear aqueous solutions (pH 7.4).		
Storage	Stored at -20°C. Avoid repeated freeze/thaw cycles.		
Application	In vitro transcription RNA amplification miRNA and siRNA synthesis		
Notes	Mix thoroughly before use and dilute with nuclease-free water if necessary.		

For Research Use Only.