

## PRODUCT INFORMATION

**GTP Solution (100 mM)**

v. 230101

<b>Catalog number</b>	C15025-1ML		
<b>Set package</b>	Cat.	Name	Amount
	C15025-1ML	GTP (100 mM)	1 mL
<b>Description</b>	Croyez Guanosine-5'-triphosphate (GTP) is supplied as a 100 mM aqueous solution at pH 7.5.		
<b>General characteristics</b>	GTP (Guanosine-5'-triphosphate, Sodium salt)		
<b>Purity</b>	≥ 98% purity of each component confirmed by HPLC.		
<b>Form</b>	Clear aqueous solutions (pH 7.4).		
<b>Storage</b>	Stored at -20°C. Avoid repeated freeze/thaw cycles.		
<b>Application</b>	In vitro transcription RNA amplification miRNA and siRNA synthesis		
<b>Notes</b>	Mix thoroughly before use and dilute with nuclease-free water if necessary.		

*For Research Use Only.*