

PRODUCT INFORMATION

Monoclonal Anti-Bet v 1 Chimeric Human IgE, clone 3B2

v. 230301

| | |
|-----------------------|---|
| Catalog number | C14005-100UG |
| Package | 100 µg |
| Product type | Monoclonal, clone 3B2 |
| Concentration | 1 mg/mL |
| Immunogen | Recombinant Bet v 1 fusion proteins |
| Isotype | Human ε, κ |
| Production | Mouse ascites |
| Preparation | KappaSelect resin purified, 0.22 µm filtered |
| Storage buffer | Phosphate-buffered saline containing 0.09 % sodium azide, pH 7.2 |
| Purity | > 95% by SDS-PAGE analysis |
| Storage | Store at -20°C. Avoid repeated freeze/thaw cycles. Suggestion: Aliquote the antibody into several vials. |
| Note | Reactivity with natural Bet v 1 proteins was tested. For research use only. Not for human in vivo or therapeutic use. |

For Research Use Only.