

PRODUCT INFORMATION

**Monoclonal Anti-Dog Serum Albumin Chimeric Human IgE,
clone 16E3**

v. 230301

Catalog number	C14011-100UG
Package	100 µg
Product type	Monoclonal, clone 16E3
Concentration	1 mg/mL
Immunogen	Dog serum albumin (Can f 3)
Isotype	Human ε, κ
Production	Mouse ascites
Preparation	KappaSelect resin purified, 0.22 µm filtered
Storage buffer	Phosphate-buffered saline containing 0.09 % sodium azide, pH 7.2
Purity	> 95% by SDS-PAGE analysis
Storage	Store at -20°C. Avoid repeated freeze/thaw cycles. Suggestion: Aliquote the antibody into several vials.
Note	Reactivity with cat serum albumin (Fel d 2) was tested. For research use only. Not for human in vivo or therapeutic use.

For Research Use Only.