

## PRODUCT INFORMATION

**Monoclonal Anti-Gad m 1 Chimeric Human IgE, clone 1F6**

v. 230301

<b>Catalog number</b>	C14017-100UG
<b>Package</b>	100 µg
<b>Product type</b>	Monoclonal, clone 1F6
<b>Concentration</b>	1 mg/mL
<b>Immunogen</b>	Pacific cod parvalbumin (Gad m 1)
<b>Isotype</b>	Human ε, κ
<b>Production</b>	Mouse ascites
<b>Preparation</b>	KappaSelect resin purified, 0.22 µm filtered
<b>Storage buffer</b>	Phosphate-buffered saline containing 0.09 % sodium azide, pH 7.2
<b>Purity</b>	> 95% by SDS-PAGE analysis
<b>Storage</b>	Store at -20°C. Avoid repeated freeze/thaw cycles. Suggestion: Aliquote the antibody into several vials.
<b>Note</b>	For research use only. Not for human in vivo or therapeutic use.

*For Research Use Only.*