


PRODUCT INFORMATION

Recombinant Der p 23 protein

v. 230301

| | |
|-----------------------|---|
| Catalog number | C14034-100UG |
| Package | 100 µg |
| Source | <i>Escherichia coli</i> |
| Sequence | ANDNDDPTTTVHPTTTEQPDDKFECPSRFGYFADPKDPHKFYICSNWEAVHK DCPGNTRWNEDEETCT with polyhistidine tag at the N-terminus |
| Concentration | 0.5 mg/mL |
| Purity | >90% as determined by SDS-PAGE. Ni-NTA chromatography |
| Form | Liquid |
| Storage buffer | Phosphate buffered saline containing 0.01% sarkosyl, pH 7.4. |
| Storage | Store at -20°C |
| Note | This product should be used and stored under sterile condition. Avoid repeated free-thaw cycles. |



SDS-PAGE analysis of Recombinant Der p 23 protein

For Research Use Only.