


PRODUCT INFORMATION

Bos d 8 protein

v. 230301

| | |
|-----------------------|---|
| Catalog number | C14050-100UG |
| Package | 100 µg |
| Concentration | 0.5 mg/mL |
| Purity | >85% as determined by SDS-PAGE. |
| Form | Liquid |
| Storage buffer | Phosphate buffered saline, pH 7.4. |
| Storage | Store at -20°C |
| Note | This product should be used and stored under sterile condition. Avoid repeated free-thaw cycles. |



SDS-PAGE analysis of Bos d 8 protein

For Research Use Only.