

## PRODUCT INFORMATION


**Gal d 2 protein**

v. 230301

---

<b>Catalog number</b>	C14047-100UG
<b>Package</b>	100 µg
<b>Concentration</b>	0.5 mg/mL
<b>Purity</b>	>95% as determined by SDS-PAGE.
<b>Form</b>	Liquid
<b>Storage buffer</b>	Phosphate buffered saline, pH 7.4.
<b>Storage</b>	Store at -20°C
<b>Note</b>	This product should be used and stored under sterile condition. Avoid repeated free-thaw cycles.

---



SDS-PAGE analysis of Gal d 2 protein

*For Research Use Only.*