

## PRODUCT INFORMATION

**Goat anti-Rabbit IgG**

v. 230201

<b>Catalog number</b>	C04041-1MG
<b>Package</b>	1 mg
<b>Description</b>	Croyez Goat anti-Rabbit IgG is purified from goat antiserum by immunoaffinity chromatography. It reacts with rabbit IgG heavy and light chains.
<b>Product type</b>	Polyclonal
<b>Concentration</b>	Lot dependent
<b>Immunogen</b>	Purified rabbit IgG
<b>Host</b>	Goat
<b>Reactivity</b>	Rabbit (Minor cross-reactivity to human in higher concentration)
<b>Conjugation</b>	N/A
<b>Isotype</b>	N/A
<b>Purity</b>	Immunogen affinity chromatography.
<b>Form</b>	Liquid
<b>Storage buffer</b>	Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.
<b>Storage</b>	Store at 4°C for at least 6 months. Store at -20°C for long-term preservation. Avoid repeated freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 4°C.
<b>Application</b>	N/A
<b>Application Note</b>	N/A
<b>Manual</b>	N/A

*For Research Use Only.*