

PRODUCT INFORMATION

Goat anti-Human IgG, Fc Fragment

v. 230201

Catalog number	C04043-1MG
Package	1 mg
Description	Croyez Goat anti-Human IgG, Fc Fragment is purified from goat antiserum by immunoaffinity chromatography. Minimal cross-reactivity with other isotype or species.
Product type	Polyclonal
Concentration	1 mg/mL
Immunogen	Gamma immunoglobulin Fc region
Host	Goat
Reactivity	Human
Conjugation	N/A
Isotype	N/A
Purity	Immunogen affinity chromatography.
Form	Liquid
Storage buffer	Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.
Storage	Store at 4°C for at least two weeks. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Application	N/A
Application Note	N/A
Manual	N/A

For Research Use Only.