

PRODUCT INFORMATION

Normal Rabbit IgG

v. 230201

Catalog number	C06004-5MG / C06004-10MG
Package	5 mg / 10 mg
Description	Rabbit IgG is purified and isolated from pooled normal rabbit serum by protein A. The purified normal rabbit IgG may be used as an isotype control in various immunoassays, such as ELISA, WB, IP, FACS, ID and IE.
Product type	Polyclonal
Concentration	Lot-dependent, please refer to product label.
Immunogen	N/A
Host	Rabbit
Reactivity	N/A
Conjugation	N/A
Isotype	N/A
Purity	>95% (SDS-PAGE)
Form	Liquid
Storage buffer	Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.
Storage	Store at 4°C for at least 6 months. Store at -20°C for long-term preservation. Avoid repeated freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 4°C.
Application	ELISA, WB, IP, FACS, ID, IE
Application Note	N/A
Manual	N/A

For Research Use Only.