

PRODUCT INFORMATION

Mouse IgG2a Kappa Isotype Control Antibody

v. 230201

Catalog number	C06007-bulk / C06007-1MG
Package	Customized package / 1 mg
Description	Mouse IgG2a Kappa isotype control antibody in used to estimate non-specific binding of target primary antibodies.
Product type	Monoclonal
Concentration	1 mg/mL
Immunogen	N/A
Host	Mouse
Reactivity	N/A
Conjugation	N/A
Isotype	IgG2a kappa
Purity	>95% (SDS-PAGE)
Form	Liquid
Storage buffer	Phosphate Buffered Saline containing 0.01% thimerosal, pH 7.4.
Storage	Store at 4°C for at least 6 months. Store at -20°C for long-term preservation. Avoid repeated freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 4°C.
Application	ELISA, WB, IP, FACS, ID, IE
Application Note	N/A
Manual	N/A

For Research Use Only.