

PRODUCT INFORMATION

Mouse anti-cMyc tag mAb sepharose

v. 230201

Catalog number	C07001-250UL / C07001-2ML / C07001-5ML
Package	0.25 mL / 2 mL / 5 mL
Description	Mouse anti-cMyc tag monoclonal antibody were conjugated to NHS-sepharose beads. Croyez Mouse anti-cMyc tag mAb sepharose has 6-10 mg mAb conjugated per mL of settled sepharose.
Product type	Monoclonal
Concentration	N/A
Immunogen	Peptide: EQKLISEEDL (cMyc)
Host	Mouse
Conjugation	Sepharose
Isotype	IgG1 kappa
Purity	N/A
Form	50% slurry
Storage buffer	Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.
Storage	Store at 4°C.
Application	IP
Application Note	N/A
Manual	The binding and elution capacity of 1 mL sepharose are commonly more than 1 mg of cMyc fusion protein. Optimal dilutions/concentrations should be determined by the end user..

For Research Use Only.