

PRODUCT INFORMATION

Goat anti-Mouse IgG sepharose

v. 230201

Catalog number	C07012-250UL / C07012-2ML / C07012-5ML
Package	0.25 mL / 2 mL / 5 mL
Description	Affinity purified antibodies were conjugated to NHS-sepharose beads. Croyez Goat anti-Mouse IgG sepharose has 6-10 mg pAb conjugated per mL of settled sepharose.
Product type	Polyclonal
Concentration	N/A
Immunogen	Purified mouse IgG
Host	Goat
Conjugation	Sepharose
Isotype	N/A
Purity	N/A
Form	50% slurry
Storage buffer	Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.
Storage	Store at 4°C.
Application	IP
Application Note	N/A
Manual	The binding and elution capacity of 1 mL sepharose are commonly more than 1 mg of IgG protein. Optimal dilutions/concentrations should be determined by the end user.

For Research Use Only.