

PRODUCT INFORMATION

Goat anti-Rabbit IgG sepharose

v. 230201

Catalog number	C07013-250UL / C07013-2ML / C07013-5ML
Package	0.25 mL / 2 mL / 5 mL
Description	Affinity purified antibodies were conjugated to NHS-sepharose. Croyez Goat anti-rabbit IgG sepharose has 6-10 mg antibodies conjugated per mL of settled sepharose.
Product type	Polyclonal
Concentration	N/A
Immunogen	Purified rabbit IgG
Host	Goat
Conjugation	Sepharose
Isotype	N/A
Purity	N/A
Form	50% slurry
Storage buffer	Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.
Storage	Store at 4°C.
Application	IP
Application Note	N/A
Manual	The binding and elution capacity of 1 mL sepharose are commonly more than 1 mg of IgG protein. Optimal dilutions/concentrations should be determined by the end user.

For Research Use Only.