

PRODUCT INFORMATION

Goat anti-Rabbit IgG (H+L)-HRP

v. 230601

Catalog number	C04003-2ML / C04003-100ML
Package	2x1 mL / 100 mL
Description	Croyez Goat anti-Rabbit IgG (H+L)-HRP is conjugated with high-activity HRP. In direct ELISA, less than 1% cross-reactivity with goat, mouse, and rat IgG is observed.
Product type	Polyclonal
Immunogen	Purified rabbit IgG
Host	Goat
Reactivity	Rabbit (Minor cross-reactivity to human in higher concentration)
Conjugation	HRP (RZ>3)
Isotype	N/A
Purity	Immunogen affinity chromatography.
Form	Liquid
Storage buffer	Phosphate Buffered Saline containing 1% BSA, 0.01% thimerosal, and 50% glycerol, pH 7.4.
Storage	Store at 4°C for at least 6 months. Store at -20°C for long-term preservation. Suggestion: Divide antibody into several vials. Keep only vials for usage at 4°C.
Application	ELISA, IHC, WB
Application Note	N/A
Manual	1:500-1:2500 for IHC 1:5000-1:25000 for ELISA and WB Working dilution for specific application should be determined by the investigator to obtain the best conditions.

For Research Use Only.