

## PRODUCT INFORMATION

**Mouse anti-M13 mAb-HRP, clone 8C1**

v. 230601

<b>Catalog number</b>	C10013-bulk / C10013-100UL
<b>Package</b>	Customized package / 100 µL
<b>Description</b>	Mouse Monoclonal antibody 8C1 recognizes bacteriophage coat protein g8p.
<b>Product type</b>	Monoclonal, clone 8C1
<b>Immunogen</b>	fd phages
<b>Host</b>	Mouse
<b>Reactivity</b>	Filamentous phages F1/fd/M13
<b>Conjugation</b>	HRP
<b>Isotype</b>	IgG <sub>2b</sub>
<b>Purification</b>	Protein A affinity chromatography
<b>Form</b>	Liquid
<b>Storage buffer</b>	HRP stabilizer
<b>Storage</b>	This product is stable after storage at: <ul style="list-style-type: none"><li>• 2-8°C for 6 months under sterile conditions from date of receipt.</li></ul> Briefly centrifuge vials before opening and protect from light.
<b>Application</b>	ELISA, Phage display
<b>Application Note</b>	N/A
<b>Manual</b>	1:250-1:1000 for ELISA Working dilution for specific application should be determined by the investigator to obtain the best conditions.

*For Research Use Only.*