

Mouse anti-M13 mAb-HRP (Antibody Set)

v. 230601

Catalog number	C10023-K01
Set package	Mouse anti-M13 mAb-HRP, clone 8-6 20 µL Mouse anti-M13 mAb-HRP, clone 11-4 20 µL Mouse anti-M13 mAb-HRP, clone 22-2 20 µL Mouse anti-M13 mAb-HRP, clone 24-3 20 µL Mouse anti-M13 mAb-HRP, clone 28-9 20 µL
Description	Mouse Monoclonal antibody recognizes bacteriophage coat protein g8p.
Product type	Monoclonal
Immunogen	fd phages
Host	Mouse
Reactivity	Filamentous phages F1/fd/M13
Conjugation	HRP
Purification	Protein A affinity chromatography
Form	Liquid
Storage buffer	Phosphate Buffered Saline containing 1% BSA, 0.03% procline 300, and 10% glycerol, pH 7.4.
Storage	This product is stable after storage at: <ul style="list-style-type: none">• 2-8°C for 6 months under sterile conditions from date of receipt. Briefly centrifuge vials before opening and protect from light.
Application	ELISA, Phage display
Application Note	N/A
Manual	1:1000-1:5000 for ELISA Working dilution for specific application should be determined by the investigator to obtain the best conditions.

For Research Use Only.