

PRODUCT INFORMATION

Goat anti-Rabbit IgG-iFluor 594

v. 230601

Catalog number	C04031-1ML / C04031-2ML
Package	1 mL / 2 mL
Description	Croyez Goat anti-Rabbit IgG-iFluor 594 is conjugated with high-activity fluorescence dye.
Product type	Polyclonal
Concentration	0.5 mg/mL
Immunogen	Purified rabbit IgG
Host	Goat
Reactivity	Rabbit (Minor cross-reactivity to human in higher concentration)
Conjugation	iFluor 594
Isotype	N/A
Purity	Immunogen affinity chromatography.
Form	Liquid
Storage buffer	Phosphate Buffered Saline containing 1% BSA, 0.01% thimerosal, pH 7.4.
Storage	<p>This product is stable after storage at:</p> <ul style="list-style-type: none"> • 2-8°C for 6 months under sterile conditions from date of receipt. <p>Do not freeze and protect from light.</p>
Application	IFA, FACS
Application Note	N/A
Manual	<p>1:200-1:1000 for IFA 1:200-1:1000 for FACS Working dilution for specific application should be determined by the investigator to obtain the best conditions.</p>

For Research Use Only.