

PRODUCT INFORMATION

Donkey anti-Goat IgG (H+L)-TAMRA

v. 230601

Catalog number	C04016-1ML / C04016-2ML
Package	1 mL / 2 mL
Description	Croyez Donkey anti-Goat IgG (H+L)-TAMRA is conjugated with high-activity fluorescence dye.
Product type	Polyclonal
Concentration	0.5 mg/mL
Immunogen	Purified goat IgG
Host	Donkey
Reactivity	Goat (Minor cross-reactivity to human in higher concentration)
Conjugation	5(6)-TAMRA
Isotype	N/A
Purity	Immunogen affinity chromatography.
Form	Liquid
Storage buffer	Phosphate Buffered Saline containing 1% BSA, 0.01% thimerosal, pH 7.4.
Storage	This product is stable after storage at: <ul style="list-style-type: none">• 2-8°C for 6 months under sterile conditions from date of receipt. Do not freeze and protect from light.
Application	IFA, FACS
Application Note	N/A
Manual	1:200-1:1000 for IFA 1:200-1:1000 for FACS Working dilution for specific application should be determined by the investigator to obtain the best conditions.

For Research Use Only.