

PRODUCT INFORMATION

UTP Solution (100 mM)

v. 240101

Catalog number	C15023-1ML		
Set package	Cat.	Name	Amount
	C15023-1ML	UTP (100 mM)	1 mL
Description	Croyez Uridine-5'-triphosphate (UTP) is supplied as a 100 mM aqueous solution at pH 7.5.		
General characteristics	UTP (Uridine-5'-triphosphate, Sodium salt)		
Purity	≥ 98% purity of each component confirmed by HPLC.		
Form	Clear aqueous solutions (pH 7.4).		
Storage	-20°C or -80°C for 12 months under sterile conditions from date of receipt.		
Application	In vitro transcription RNA amplification miRNA and siRNA synthesis		
Notes	Mix thoroughly before use and dilute with nuclease-free water if necessary.		

For Research Use Only.