


PRODUCT INFORMATION

Recombinant Bla g 5 protein

v. 231201

Catalog number	C14055-100UG
Package	100 µg
Source	<i>Escherichia coli</i>
Sequence	APSYKLTYPVKALGEPFRLLSYGEKDFEDYRFQEGDWPNLKPSMPFGKTPV LEIDGKQTHQSVAISRYLGKQFGLSGKDDWENLEIDMIVDTISDFRAAIANYHYD ADENSKQKKWDPLKKETIPYYTKKFDEVVKANGGYLAAGKLTWADFYFVAILDY LNHMAKEDLVANQPNLKALREKVLGLPAIKAWVAKRPPTDL with polyhistidine tag at the N-terminus
Concentration	0.5 mg/mL
Purity	>90% as determined by SDS-PAGE. Ni-NTA chromatography
Form	Liquid
Storage buffer	Phosphate buffered saline, pH 7.4.
Storage	Store at -20°C
Note	This product should be used and stored under sterile condition. Avoid repeated free-thaw cycles.



SDS-PAGE analysis of Recombinant Bla g 5 protein

For Research Use Only.