

PRODUCT INFORMATION

Mouse Anti- Human IgE- HRP

v. 240401

Catalog number	C04051-1ML
Package	1 mL
Description	Mouse Anti- Human IgE-HRP is conjugated with high-activity horseradish peroxidase (HRP).
Clonality	Monoclonal
Clone Name	Clone 102-1
Host	Mouse
Reactivity	Human IgE
Conjugation	HRP
Form	Liquid
Storage buffer	Phosphate Buffered Saline containing 1% BSA, 0.03% ProClin 300 and 10% glycerol, pH 7.4
Stability & Storage	<p>This product is stable after storage at:</p> <ul style="list-style-type: none">• 2-8°C for 2 weeks under sterile conditions from date of receipt.• -20°C or -80°C for 12 months under sterile conditions from date of receipt. <p>Avoid repeated freeze/thaw cycles and protect from light.</p>
Application	ELISA
Product Note	<ol style="list-style-type: none">1. Recommended dilution factor: ELISA: 2000-10000 <p>Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.</p> <ol style="list-style-type: none">2. Divide antibody into several vials. Keep only vials for usage at 2-8°C.

For Research Use Only.