

## PRODUCT INFORMATION

**Monoclonal Anti-Gal d 2 Chimeric Human IgE, clone 9G5**

v. 250401

|                       |   |
|-----------------------|---|
| <b>Catalog number</b> | C14019-100UG  |
| <b>Package</b>        | 100 µg  |
| <b>Product type</b>   | Monoclonal, clone 9G5   |
| <b>Concentration</b>  | 1 mg/mL   |
| <b>Immunogen</b>      | Ovalbumin   |
| <b>Isotype</b>        | Human ε, κ  |
| <b>Production</b>     | Mouse ascites   |
| <b>Preparation</b>    | KappaSelect resin purified, 0.22 µm filtered  |
| <b>Storage buffer</b> | Phosphate-buffered saline containing 0.09 % sodium azide, pH 7.2  |
| <b>Purity</b>         | > 95% by SDS-PAGE analysis  |
| <b>Storage</b>        | Store at -20°C. Avoid repeated freeze/thaw cycles.<br>Suggestion: Aliquote the antibody into several vials. |
| <b>Note</b>           | For research use only. Not for human in vivo or therapeutic use.  |

*For Research Use Only.*