

PRODUCT INFORMATION

Anti-Influenza B virus NP Antibody [clone IB06]

v. 251001

Catalog number	C10055-bulk / C10055-100UG
Package	Customized package / 100 µg
Description	Anti-Influenza B virus NP Antibody [clone IB06] only recognizes nucleocapsid protein (NP) of influenza B viruses but not influenza A viruses. Influenza B viruses are RNA viruses and divided into two genetic lineages (B/Yamagata and B/Victoria). Influenza B viruses are not classified by subtype like influenza A viruses.
Clonality	Recombinant Mouse IgG
Reactivity	Influenza B virus
Conjugation	Unconjugated
Isotype	IgG _{2a}
Concentration	1 mg/mL
Form	Liquid
Storage buffer	Phosphate Buffered Saline, pH 7.4.
Stability & Storage	This product is stable after storage at: • 2-8°C for 2 weeks under sterile conditions from date of receipt. • -20°C or -80°C for 12 months under sterile conditions from date of receipt. Avoid repeated freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.
Application	ELISA, CLIA, LFIA

For Research Use Only.