

Achieve High-Efficiency HRP Conjugation in Minutes with the Croyez HRP Conjugation Kit



Ideal for small-scale antibody or protein labelling (10 µg–1 mg) for ELISA, Western blotting and IHC

Croyez simplifies the complex chemical conjugation process into three easy steps — mix, incubate, and complete — enabling faster, more stable, and highly efficient HRP labeling.

The kit achieves high-activity HRP conjugation in a short time, ensuring both enzyme activity and labeling stability, and is ideal for ELISA, Western blot, and immunohistochemistry (IHC) applications.

Rapid & easy:

Full kit includes high-activity HRP and all reagents for conjugation.

Broad application:

Compatible with ELISA, Western blot, IHC.


Optimised for small quantities:

Designed for 10 µg to 1 mg of antibody or protein.

Flexible buffer tolerance: Up to 20 mM Tris, 10 % glycerol, 0.02 % NaN₃ have minimal effect.

Stable and user-friendly storage: –20 °C or –80 °C, 12 months under sterile conditions.

How it works



Cat#	Product	Package	Cat#	Product	Package
C08001-10	HRP Conjugation Kit	10 µg/set	C04005-1ML	Rabbit anti-Mouse IgG (H+L)-HRP	1 mL
C08001-100	HRP Conjugation Kit	100 µg/set	C04005-2ML	Rabbit anti-Mouse IgG (H+L)-HRP	2 mL
C08001-1000	HRP Conjugation Kit	1 mg/set	C04006-1ML	Donkey anti-Goat IgG (H+L)-HRP	1 mL
C04051-1ML	Mouse Anti-Human IgE-HRP	1 mL	C04006-2ML	Donkey anti-Goat IgG (H+L)-HRP	2 mL
C10013-100UL	Mouse anti-M13 mAb-HRP, clone 8C1	100 µL	C04007-1ML	Sheep anti-Chicken IgY-HRP	1 mL
C10014-100UL	Mouse anti-M13 mAb-HRP, clone 8-6	100 µL	C04007-2ML	Sheep anti-Chicken IgY-HRP	2 mL
C04001-2ML	Goat anti-Mouse IgG (H+L)-HRP	2x1 mL	C04009-1ML	AccuBlot™ anti-Mouse IgG-HRP	1 mL
C04002-1ML	Goat anti-Rat IgG (H+L)-HRP	1 mL	C04009-2ML	AccuBlot™ anti-Mouse IgG-HRP	2 mL
C04002-2ML	Goat anti-Rat IgG (H+L)-HRP	2 mL	C04010-1ML	AccuBlot™ anti-Rabbit IgG-HRP	1 mL
C04003-2ML	Goat anti-Rabbit IgG (H+L)-HRP	2x1 mL	C04010-2ML	AccuBlot™ anti-Rabbit IgG-HRP	2 mL
C04004-1ML	Goat anti-Human IgG (Fc Fragment)-HRP	1 mL	C08009-500UL	Streptavidin-HRP	500 µL
C04004-2ML	Goat anti-Human IgG (Fc Fragment)-HRP	2 mL			

